

**Color coded
size indicator -**
push in until flush
with handle to
activate batteries

1



2

On/off button -
click to activate
radial LED light
source

ONETRAC® LXS CONTROL GUIDE



ONETRAC LXS is a single-use suction device with **integrated, cordless radial LED light source**.

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the surgical setting by making reusable devices and fiber-optic light cables obsolete.

Illumination and aspiration in the depths of a surgical pocket can be met within seconds – no assembly, additional parts/components, or fiber-optic cables required.

To learn more visit www.obpmedical.com/lxs

ONETRAC® LXS

Single-Use Suction Device with Integrated Cordless Radial LED Light Source

SINGLE-USE CORDLESS ILLUMINATION